




Urology Instruments

Urologische Instrumente



EM-13-2022-01 19.5 cm
EM-13-2022-02 22 cm
EM-13-2022-03 23 cm



V1



PEAN
EM-13-2022-01 - 2022-03



V2

GUYON
EM-13-2023
23.5 cm

Atraumatic



EM-13-2025-02
24 cm

EM-13-2025-01
23 cm



GUYON
EM-13-2024
23 cm



GUYON
EM-13-2025-01 - 2025-03




EM-13-2025-03
23 cm



STILLE
EM-13-2026
22 cm



STILLE
EM-13-2027
24 cm



GUYON
EM-13-2028
26.5 cm



HERRICK
EM-13-2029
22.5 cm



WERTHEIM
EM-13-2030
25.5 cm

Atraumatic



WERTHEIM - CULLEN
EM-13-2031
Angle- and pressure clamp
21.5 cm



WERTHEIM
EM-13-2032
Angle- and pressure clamp
25 cm



1/2

DARTIGUES
EM-13-2033
25.5 cm



1/2

COLLIN
EM-13-2034
28.5 cm



EM-13-2035
17 cm

EM-13-2036
16 cm



EM-13-2037

EM-13-2038




HEIMING
EM-13-2035 - 2036
Çocuklar için/For children



RANDALL
EM-13-2037 - 2040
22.5 cm



Vi
EM-13-2041



Vi
EM-13-2044
24.5 cm



Vi
EM-13-2042




Vi
EM-13-2043
23 cm



Vi
RANDALL
EM-13-2041 - 2042
19 cm



RANDALL
EM-13-2043 - 2044
Extremely delicate patterns



EM-13-2045-04
57 x 25 mm



EM-13-2045-04
57 x 25 mm




EM-13-2045-05
82 x 25 mm



EM-13-2045-05
82 x 25 mm

- MILLIN
EM-13-2045-01 Complete
EM-13-2045-03 Retractor only
EM-13-2045-04 Blades 57 x 25 mm
EM-13-2045-05 Blades 82 x 25 mm
EM-13-2045-06 Center blade 105 x 45 mm


- MILLIN
EM-13-2045-02 Complete
With light carrier
EM-13-2045-03 Retractor only
EM-13-2045-04 Blades 57 x 25 mm
EM-13-2045-05 Blades 82 x 25 mm
EM-13-2045-07 Center blade 105 x 45 mm



67 mm

1/2

KOCHER
EM-13-2046
Bladder retractor
26 cm



90 mm

35 mm

1/2

LEGUEU
EM-13-2047
Bladder spatula
26 cm



1/2


1/2

MILLIN
EM-13-2048
Bladder neck spreader
28.5 cm

Atraumatic



Atraumatic



MILLIN
EM-13-2049
Lobe grasping forceps
23.5 cm



MILLIN
EM-13-2050
Lobe grasping forceps
24.5 cm



MILLIN
EM-13-2051
Ligature holding forceps
24.5 cm



1/2

YOUNG
EM-13-2052
21 cm



1/2

YOUNG
EM-13-2053
21 cm




1/2

EM-13-2054-01
Prostate retractor
36 cm



1/2

EM-13-2054-02
Urethral retractor
39 cm



EM-13-2055
Venous plexus clamp
21 cm

BARRE
EM-13-2056-01 - 2056-02
Dissecting forceps
28 cm

BARRE
EM-13-2057
Dissecting scissors
28 cm



BARRE
EM-13-2058
No.3 XL
For left hand
25 cm



BARRE
EM-13-2059-01 - 2059-02
No.3 XL
For left hand
25 cm



EM-13-2059-01 12°
EM-13-2059-02 25°



YOUNG-HRYNTSCHAK
EM-13-2060
With 1 needle F 400.00
24 cm



F 400.00
Vi



Vi
EM-13-2061-01



Vi
EM-13-2061-02



Ucu tam izoleli
Fully insulated

BARRE
EM-13-2062
Bipolar forceps
26 cm



EM-13-2063
Ø 4.8 mm
Forceps for grasping the glans
and the foreskin
12.5 cm




EM-13-2064
Fixation clamp for the
spermatic cord
14.5 cm

Atraumatic



EM-13-2065
Hook for the spermatic cord
16.5 cm




EM-13-2066
Ureter forceps
24 cm

- Foreskin is pushed over the bell-shaped protection rod with the help of two hemostatic or ligature forceps and the glans is placed into cavity on the rod.
- In this position, wing nut is tightened as close as possible for guaranteeing annular copression on the foreskin. Then foreskin is cut off with scalpel, right next to the compression ring.
- The circumcision clamp is kept in compressed position for few minutes for avoiding bleeding. In case of large vein of frenular artery is cut artrey is cut through, the ligature set should be kept ready.
- The advantage of the method described is aesthetically favorble scar after opretion.



WINKELMANN
EM-13-2067-01 Ø 12 mm
EM-13-2067-02 Ø 16 mm
EM-13-2067-03 Ø 20 mm
Circumcision clamps



STOCKMANN
EM-13-2068
Penis klempe
Penis clamp
7 cm



STRAUSS
EM-13-2069
Penis klempe
Penis clamp
13 cm



MATHIEU
EM-13-2070
Foreign body forceps
19 cm



PITHA
EM-13-2071
Foreign body forceps
22 cm



ASCH
EM-13-2072
Uterine secretion scoop, blunt
20 cm



EM-13-2073-03 - 2073-06
Uterine probe
33 cm

| | Charr. |
|---------------|--------|
| EM-13-2073-03 | 3 |
| EM-13-2073-04 | 4 |
| EM-13-2073-05 | 5 |
| EM-13-2073-06 | 6 |



1/2

| | Charr. |
|---------------|--------|
| EM-13-2074-08 | 8 |
| EM-13-2074-10 | 10 |
| EM-13-2074-12 | 12 |
| EM-13-2074-14 | 14 |
| EM-13-2074-16 | 16 |
| EM-13-2074-18 | 18 |
| EM-13-2074-20 | 20 |

EM-13-2074-08 - 2074-20
Male catheters, metal
29 cm



1/1

| | Charr. |
|---------------|--------|
| EM-13-2075-06 | 6 |
| EM-13-2075-07 | 7 |
| EM-13-2075-08 | 8 |
| EM-13-2075-09 | 9 |
| EM-13-2075-10 | 10 |
| EM-13-2075-11 | 11 |
| EM-13-2075-12 | 12 |
| EM-13-2075-13 | 13 |
| EM-13-2075-14 | 14 |
| EM-13-2075-15 | 15 |
| EM-13-2075-16 | 16 |
| EM-13-2075-18 | 18 |
| EM-13-2075-20 | 20 |

EM-13-2075-06 - 2075-20
Female catheters, metal
15.5 cm



1/4

EM-13-2076
Catheter guide
49 cm



1/3
GUYON
EM-13-2077
Catheter guide
36.5 cm



1/3
EM-13-2078
Introducer for catheters
46 cm



1/2
DITTEL
EM-13-2079-08 - 2079-30
Urethral dilating bougies
34.5 cm

| | Charr. |
|---------------|--------|
| EM-13-2079-08 | 8 |
| EM-13-2079-09 | 9 |
| EM-13-2079-10 | 10 |
| EM-13-2079-11 | 11 |
| EM-13-2079-12 | 12 |
| EM-13-2079-13 | 13 |
| EM-13-2079-14 | 14 |
| EM-13-2079-15 | 15 |
| EM-13-2079-16 | 16 |
| EM-13-2079-17 | 17 |
| EM-13-2079-18 | 18 |
| EM-13-2079-19 | 19 |
| EM-13-2079-20 | 20 |
| EM-13-2079-21 | 21 |
| EM-13-2079-22 | 22 |
| EM-13-2079-23 | 23 |
| EM-13-2079-24 | 24 |
| EM-13-2079-25 | 25 |
| EM-13-2079-26 | 26 |
| EM-13-2079-27 | 27 |
| EM-13-2079-28 | 28 |
| EM-13-2079-29 | 29 |
| EM-13-2079-30 | 30 |

| | |
|---------------|-------------------|
| EM-13-2080-24 | Ben 24/ Charr. 12 |
| EM-13-2080-25 | Ben.25 |
| EM-13-2080-26 | Ben 26/ Charr. 13 |
| EM-13-2080-27 | Ben.27 |
| EM-13-2080-28 | Ben 28/ Charr. 14 |
| EM-13-2080-29 | Ben.29 |
| EM-13-2080-30 | Ben 30/ Charr. 15 |
| EM-13-2080-31 | Ben.31 |
| EM-13-2080-32 | Ben 32/ Charr. 16 |
| EM-13-2080-33 | Ben.33 |
| EM-13-2080-34 | Ben 34/ Charr. 17 |
| EM-13-2080-35 | Ben.35 |
| EM-13-2080-36 | Ben 36/ Charr. 18 |
| EM-13-2080-37 | Ben.37 |
| EM-13-2080-38 | Ben 38/ Charr. 19 |
| EM-13-2080-39 | Ben.39 |
| EM-13-2080-40 | Ben 40/ Charr. 20 |
| EM-13-2080-41 | Ben.41 |
| EM-13-2080-42 | Ben 42/ Charr. 21 |
| EM-13-2080-43 | Ben.43 |
| EM-13-2080-44 | Ben 44/ Charr. 22 |
| EM-13-2080-45 | Ben.45 |
| EM-13-2080-46 | Ben 46/ Charr. 23 |
| EM-13-2080-47 | Ben.47 |
| EM-13-2080-48 | Ben 48/ Charr. 24 |
| EM-13-2080-49 | Ben.49 |
| EM-13-2080-50 | Ben 50/ Charr. 25 |
| EM-13-2080-51 | Ben.51 |
| EM-13-2080-52 | Ben 52/ Charr. 26 |
| EM-13-2080-53 | Ben.53 |
| EM-13-2080-54 | Ben 54/ Charr. 27 |
| EM-13-2080-55 | Ben.55 |
| EM-13-2080-56 | Ben 56/ Charr. 28 |
| EM-13-2080-57 | Ben.57 |
| EM-13-2080-58 | Ben 58/ Charr. 29 |
| EM-13-2080-59 | Ben.59 |
| EM-13-2080-60 | Ben 60/ Charr. 30 |



1/2

GUYON
EM-13-2080-24 - 2080-60
Urethral dilating bougies
25.5 cm



1/2

TURNER WARWICK
EM-13-2081
Retrograde bougie
18 cm