



EXPO-MED

Serving Health & Medical Division Since 1990

06

Parodontal Periodontal



EM-3377-1 For anterior incisors and canines, slightly angled
GRA 1/2



EM-3377-3 For anterior and premolars, slightly angled
GRA 3/4



EM-3377-5 For anterior incisors and canines, angled
GRA 5/6



EM-3377-7 For premolars/molars, buccal/lingual, strongly angled
GRA 7/8



EM-3377-9 For premolars/molars, buccal/lingual, very strongly angled
GRA 9/10



EM-3377-11 For all mesial surfaces of premolars/molars
GRA 11/12






EM-3377-13 For all distal surfaces of premolars/molars
GRA 13/14



EM-3377-15 For premolars/molars, all mesial surfaces, angulation like 13/14, but opposite blades
GRA 15/16



EM-3377-17 For premolars/molars, for optimal access to distal surfaces and deep pockets
GRA 17/18



EM-3377-101 Deep scaling for anteriors and premolars, angled
GRA 1/2M5



EM-3377-103 Deep scaling for anteriors and premolars, angled
GRA 3/4M5



EM-3377-105 Deep scaling for anteriors and premolars, sharply angled
GRA 5/6M5



EM-3377-107 Deep scaling for premolars and molars, buccal/lingual, angled
GRA 7/8M5



EM-3377-109 Deep scaling for premolars and molars, buccal/lingual, sharply angled
GRA 9/10M5



EM-3377-111 Deep scaling for premolars and molars, mesial, double angled
GRA 11/12M5



EM-3377-113 Deep scaling for premolars and molars, distal, double angled
GRA 13/14M5



EM-3303-01TI L 1/2 LANGER L 1/2, universal curette, titanium tips, for mandibular posterior




EM-3303-03TI L 3/4 LANGER L 3/4, universal curette, titanium tips, for maxillary posterior



EM-3303-05TI L 5/6 LANGER L 5/6, universal curette, titanium tips, for maxillary and mandibular anterior



EM-3302-04TI CU 4R/4L COLUMBIA UNIVERSITY 4R/4L, universal curette, titanium tips, for maxillary and mandibular posterior



EM-3300-00 Key to exchange working tips



EM-3302-02 COLUMBIA UNIVERSITY 2R/2L, for front teeth
CU 2R/2L



EM-3302-13CU COLUMBIA UNIVERSITY 13/14, universal curette, especially for root planing
CU 13/14



EM-3302-04 COLUMBIA UNIVERSITY 4R/4L, universal curette, for maxillary and mandibular posterior
CU 4R/4L



EM-3302-13 MCCALL 13S/14S, universal scaler, for interdental spaces of front teeth and premolars, fine, sharply pointed tips
MC 13S/14S



EM-3302-17 MCCALL 17S/18S, universal scaler for molars, fine, rounded tips
MC 17S/18S



EM-3303-01 LANGER L 1/2, universal curette, for mandibular posterior
L 1/2




EM-3303-03 LANGER L 3/4, universal curette, for maxillary posterior
L 3/4



EM-3303-05 LANGER L 5/6, universal curette, for maxillary and mandibular anterior
L 5/6



EM-3302-23A Universal curette M 23A, for upper and lower jaw teeth and supragingival calculus, slim
M23A



EM-3303-100

Prophylaxe Kit Universal, consisting of:
 Universal curettes Langer 1/2, 3/4, 5/6, M23A,
 Scaler 204S and tray 2599-05



EM-3377-91 L 1/2 LANGER L 1/2, universal curette, for mandibular posterior



EM-3377-93 L 3/4 LANGER L 3/4, universal curette, for maxillary posterior



EM-3377-95 L 5/6 LANGER L 5/6, universal curette, for maxillary and mandibular anterior



EM-3305-06 HYGIENIST sickle scaler with opposed tips, for anterior and premolars, mesial/distal
H 6/7



EM-3305-533 HYGIENIST/JACQUETTE Scaler 5/33, for anterior and premolars, especially for heavy calculus
H5/33



EM-3305-02 KAPLAN sickle scaler with large facial surfaces matched to tooth contours for removing heavy calculus, plastic, cement, or adhesive residues from all interdental spaces
CI 2/3



EM-3305-35 Scaler with small facial surfaces matched to tooth contours, for all interdental spaces
U135



EM-3306-02 TAYLOR medium-sized scaler with facial surfaces matched to tooth contours, for all interdental spaces
T2/3



EM-3305-04S Sickle scaler for removing dental plaque from interdental spaces in the molar area
204S



EM-3305-04SD Sickle scaler for removing dental plaque from narrow interdental spaces in the premolar and frontal area.
204SD



EM-3305-23 Universal sickle scaler for removing calculus from interdental spaces
M23



EM-3305-GX Universal sickle curette (scaler), rounded tips
GXC4



EM-3253 CUMINE



EM-3254 MITCHELL Osteotrimmer



EM-3274-204
2045 Sickle scaler for removing dental plaque from interdental spaces in the molar area



EM-3274-6 HYGIENIST sickle scaler with opposed tips, for anterior and premolars, mesial/distal
H6/H7



EM-3275-23 TAYLOR medium-sized scaler with facial surfaces matched to tooth contours, for all interdental spaces
T2/3



EM-3276 CRANE-KAPLAN 2/3
CK6



EM-3380-131 MCCALL 13S/14S, universal scaler, for interdental spaces of front teeth and premolars, fine, sharply pointed tips
13S/14S



EM-3274-533 HYGIENIST/JACQUETTE Scaler 5/33, for anterior and premolars, especially for heavy calculus
H5/33




EM-3376-5 GOLDMANN-FOX hoe for removing heavy calculus deposits, buccal/lingual
GF5



EM-3376-6 GOLDMANN-FOX hoe for mesial and distal surfaces of the posterior teeth
GF6

EM-3308-100

Furcation set for removal of calculus and infected soft tissue within furcations.
Tungsten coated working tips.



EM-3308-01 Furcation Curette buccal/lingual 0.9 mm



EM-3308-02 Furcation Curette buccal/lingual 1.3 mm



EM-3308-03 Furcation Curette mesial/distal 0.9 mm



EM-3308-04 Furcation Curette mesial/distal 1.3 mm



EM-3308-05 Furcation Scaler labial-lingual



EM-3395
15K/16K KIRKLAND particularly used in retromolar areas



EM-3396
GF7 GOLDMANN-FOX for long, bevel incisions, their teardrop-shaped tips suite also in retromolar areas



EM-3397
GF8 GOLDMANN-FOX spear-shaped, knife for interdental areas



EM-3398
GF11 GOLDMANN-FOX for molar areas



EM-3399
1/2 ORBAN for preparation in molar areas



EM-3232
3/7 HIRSCHFELD FIG. 3/7, buccal/lingual



EM-3232-2
5/11 HIRSCHFELD FIG. 5/11, mesial/distal



EM-3232-6
35/4 SUGARMAN FIG. 35/4, buccal/lingual



EM-3232-7
15/25 SUGARMAN FIG. 15/25, mesial/distal




EM-3233-9
9/10 SCHLUGER FIG. 9/10, curved



EM-3239-10
COHEN
10,0 cm

EM-3239-13
ILIC
14,0 cm



EM-4085-00TI

Composite set, consisting of: 4085-01Ti, 4085-02Ti, 4085-03Ti, 4085-04Ti, tray 2599-05



EM-4085-01TI Composite spatula FIG. 1, fine, 1,8 / 1,5 mm, exchangeable



EM-4085-02TI Composite spatula FIG. 2, 2,9 / 2,7 mm, exchangeable



EM-4085-03TI Composite pear shape plugger, Ø 1,5 / 2,0 mm, exchangeable



EM-4085-04TI Composite occlusal former / burnisher, Ø 1,8 / 1,8 mm, exchangeable, one tip angulated especially for distal



EM-4085-05Ti Composite ball plugger, fine, Ø 0,9 / 1,3 mm, exchangeable



EM-4085-06Ti Composite ball plugger, Ø 2,1 / 1,3 mm flat, exchangeable




EM-4085-07Ti Composite plugger, Ø 1,6 / 2,5 mm, exchangeable



EM-4085-08Ti Composite spatula, 1,7 / 1,7mm, exchangeable, especially interproximal mesial / distal



EM-4086-10 Composite spatula, fine, 1,6 / 1,6 mm, rigid



EM-4086-110 Composite spatula, 2,0 / 2,0 mm, rigid



EM-4086-1 HEIDEMANN FIG. 0 composite spatula, 2,0 / 2,0 mm, flexible



EM-4086-2 HEIDEMANN FIG. 1 composite spatula, 2,5 / 2,5 mm, flexible



EM-4086-3 HEIDEMANN FIG. 2 composite spatula, 3,0 / 3,0 mm, flexible