



EXPO-MED

Serving Health & Medical Division Since 1990

09

Abdrucklöffel

Impression



9. Abdrucklöffel · Impression Trays

Regular Pattern Rimlock

INOX
STAINLESS STEEL

U1/L1  XL

U2/L2  L

U3/L3  M

U4/L4  S

U5/L5  XS



EM-8001-U1

EM-8001-U2

EM-8001-U3



EM-8001-U4

EM-8001-U5



		A mm.	B mm.
EM-8001-U1	XL	81	63
EM-8001-U2	L	74	60
EM-8001-U3	M	71	55
EM-8001-U4	S	67	53
EM-8001-U5	XS	65	51

		A mm.	B mm.
EM-8001-L1	XL	80	63
EM-8001-L2	L	78	59
EM-8001-L3	M	76	57
EM-8001-L4	S	72	57
EM-8001-L5	XS	69	50



EM-8001-L1

EM-8001-L2

EM-8001-L3

EM-8001-L4

EM-8001-L5



9. Abdrucklöffel · Impression Trays

Regular Pattern Rimlock

INOX
STAINLESS STEEL

U1/L1 ○ XXS

U2/L2 ● XS

U3/L3 ● S

U4/L4 ● M

U5/L5 ● L

U6/L6 ● XL



EM-8002-U6

EM-8002-U5

EM-8002-U4

EM-8002-U3



EM-8002-U2

EM-8002-U1



		A	B
		mm.	mm.
EM-8002-U1	XXS	61	47
EM-8002-U2	XS	64	51
EM-8002-U3	S	66	53
EM-8002-U4	M	69	57
EM-8002-U5	L	71	60
EM-8002-U6	XL	73	62

		A	B
		mm.	mm.
EM-8002-L1	XXS	55	47
EM-8002-L2	XS	60	49
EM-8002-L3	S	72	54
EM-8002-L4	M	74	57
EM-8002-L5	L	76	59
EM-8002-L6	XL	78	61



XL

L

M

S

XS

XXS



EM-8002-L6

EM-8002-L5

EM-8002-L4

EM-8002-L3

EM-8002-L2

EM-8002-L1



EM-8002-L2D

EM-8002-L4D

EM-8002-L6D



EM-8002-U1P

EM-8002-L1P



EM-8002-L2S

EM-8002-L4S

EM-8002-L6S



9. Abdrucklöffel · Impression Trays

Regular Pattern Without Rimlock




INOX
STAINLESS STEEL

- U1/L1 ○ XXS
- U2/L2 ● XS
- U3/L3 ● S
- U4/L4 ● M
- U5/L5 ● L
- U6/L6 ● XL






EM-8003-U6 EM-8003-U5 EM-8003-U4 EM-8003-U3 EM-8003-U2 EM-8003-U1



		A mm.	B mm.
EM-8003-U1	XXS	61	47
EM-8003-U2	XS	64	51
EM-8003-U3	S	66	53
EM-8003-U4	M	69	57
EM-8003-U5	L	71	60
EM-8003-U6	XL	73	62

		A mm.	B mm.
EM-8003-L1	XXS	55	47
EM-8003-L2	XS	60	49
EM-8003-L3	S	72	54
EM-8003-L4	M	74	57
EM-8003-L5	L	76	59
EM-8003-L6	XL	78	61





9. Abdruckklöffel · Impression Trays

Regular Pattern Perforated Rimlock



INOX
STAINLESS STEEL

U1/L1 ○ XXS

U2/L2 ● XS

U3/L3 ● S

U4/L4 ● M

U5/L5 ● L

U6/L6 ● XL





9. Abdruckklöffel · Impression Trays

Regular Pattern Perforated Rimlock

EM-8004-U6


EM-8004-U5

EM-8004-U4

EM-8004-U3




EM-8004-U2

EM-8004-U1



		A mm.	B mm.
EM-8004-U1	XXS	61	47
EM-8004-U2	XS	64	51
EM-8004-U3	S	66	53
EM-8004-U4	M	69	57
EM-8004-U5	L	71	60
EM-8004-U6	XL	73	62

		A mm.	B mm.
EM-8004-U1	XXS	55	47
EM-8004-U2	XS	60	49
EM-8004-U3	S	72	54
EM-8004-U4	M	74	57
EM-8004-U5	L	76	59
EM-8004-U6	XL	78	61





9. Abdruckklöffel · Impression Trays

Regular Pattern Perforated Without Rimlock



INOX
STAINLESS STEEL

U1/L1 ○ XXS

U2/L2 ● XS

U3/L3 ● S

U4/L4 ● M

U5/L5 ● L

U6/L6 ● XL





9. Abdruckklöffel · Impression Trays

Regular Pattern Perforated Without Rimlock

EM-8005-U6

EM-8005-U5

EM-8005-U4

EM-8005-U3

EM-8005-U2

EM-8005-U1



XL

L

M




S

XS

XXS

		A mm.	B mm.
EM-8005-U1	XXS	61	47
EM-8005-U2	XS	64	51
EM-8005-U3	S	66	53
EM-8005-U4	M	69	57
EM-8005-U5	L	71	60
EM-8005-U6	XL	73	62

		A mm.	B mm.
EM-8005-L1	XXS	55	47
EM-8005-L2	XS	60	49
EM-8005-L3	S	72	54
EM-8005-L4	M	74	57
EM-8005-L5	L	76	59
EM-8005-L6	XL	78	61



EM-8005-L6


EM-8005-L5

EM-8005-L4

EM-8005-L3

EM-8005-L2

EM-8005-L1



EM-8005-L2D

EM-8005-L4D

EM-8005-L6D

EM-8005-L2S

EM-8005-L4S

EM-8005-L6S



9. Abdrucklöffel · Impression Trays

Regular Pattern Perforated Aluminum Without Rimlock





9. Abdruckklöffel · Impression Trays

Regular Pattern Perforated Aluminum Without Rimlock

EM-8006-U6

EM-8006-U5

EM-8006-U4

EM-8006-U3


EM-8006-U2

EM-8006-U1



		A mm.	B mm.
EM-8006-U1	XXS	63	48
EM-8006-U2	XS	66	51
EM-8006-U3	S	68	54
EM-8006-U4	M	71	58
EM-8006-U5	L	74	62
EM-8006-U6	XL	76	63

		A mm.	B mm.
EM-8006-L1	XXS	56	48
EM-8006-L2	XS	59	50
EM-8006-L3	S	72	52
EM-8006-L4	M	76	57
EM-8006-L5	L	77	56
EM-8006-L6	XL	78	59



XL

L

M

S

XS

XXS



EM-8006-L6

EM-8006-L5

EM-8006-L4

EM-8006-L3

EM-8006-L2

EM-8006-L1



EM-8006-L2D

EM-8006-L4D

EM-8006-L6D

EM-8006-L2S

EM-8006-L4S

EM-8006-L6S



9. Abdruckklöffel · Impression Trays

Premium Plastic Molded Material



EM-8007-U4

EM-8007-U3

EM-8007-U2

EM-8007-U1

EM-8007-U0



XL

L



M

S

XS

		A mm.	B mm.
EM-8007-U0	XS	65	51
EM-8007-U1	S	67	53
EM-8007-U2	M	71	55
EM-8007-U3	L	74	60
EM-8007-U4	XL	81	63

		A mm.	B mm.
EM-8007-L0	XS	69	50
EM-8007-L1	S	72	57
EM-8007-L2	M	76	57
EM-8007-L3	L	78	59
EM-8007-L4	XL	80	63



XL

L

M

S

XS



EM-8007-L4

EM-8007-L3

EM-8007-L2

EM-8007-L1

EM-8007-L0



EM-8008-U6

EM-8008-U5

EM-8008-U4

EM-8008-U3

EM-8008-U2

EM-8008-U1



XL

L

M

S

XS

XXS

EM-8008-L6

EM-8008-L5

EM-8008-L4

EM-8008-L3

EM-8008-L2

EM-8008-L1



XL

L

M

S

XS

XXS



EM-8009-U7

EM-8009-U4


EM-8009-U1




EM-8009-L22

EM-8009-L21

EM-8009-L20

sizes		
	A mm.	B mm.
EM-8009-U1	L 81	67
EM-8009-U4	M 71	64
EM-8009-U7	S 70	57

sizes		
	A mm.	B mm.
EM-8009-L20	L 81	60
EM-8009-L21	M 79	57
EM-8009-L22	S 76	54



EM-8010-U7



EM-8010-U4



EM-8010-U1





EM-8010-L22





EM-8010-L21





EM-8010-L20


sizes			
			
	A	B	
	mm.	mm.	
EM-8010-U1	L	81	67
EM-8010-U4	M	71	64
EM-8010-U7	S	70	57
sizes			
			
	A	B	
	mm.	mm.	
EM-8010-L20	L	81	60
EM-8010-L21	M	79	57
EM-8010-L22	S	76	54





sizes			sizes				
	A mm.	B mm.		A mm.	B mm.		
EM-8011-U61	L	79	70	EM-8011-L66	L	86	64
EM-8011-U62	M	76	67	EM-8011-L67	NL	70	64
EM-8011-U63	NM	70	64	EM-8011-L68	M	73	60
EM-8011-U64	S	73	60	EM-8011-L69	S	70	57



sizes			sizes				
	A	B		A	B		
	mm.	mm.		mm.	mm.		
EM-8012-U61	L	79	70	EM-8012-L66	L	86	64
EM-8012-U62	M	76	67	EM-8012-L67	NL	70	64
EM-8012-U63	NM	70	64	EM-8012-L68	M	73	60
EM-8012-U64	S	73	60	EM-8012-L69	S	70	57




 <small>sizes</small>			A	B
			mm.	mm.
	EM-8013-U1	XXS	60	50
	EM-8013-U2	XS	64	54
	EM-8013-U3	S	68	58
	EM-8013-U4	M	70	60
	EM-8013-U5	L	78	64

 <small>sizes</small>			A	B
			mm.	mm.
	EM-8013-L1	XXS	60	54
	EM-8013-L2	XS	67	55
	EM-8013-L3	S	70	56
	EM-8013-L4	M	74	58
	EM-8013-L5	L	76	60
	EM-8013-L6	XL	80	62




For upper jaws
for children

	A	B
	mm.	mm.
EM-8014-OB 0	59	52
EM-8014-OB 1	66	52
EM-8014-OB 2	71	62
EM-8014-OB 3	71	64
EM-8014-OB 4	80	70



For lower jaws
for children




	A	B
	mm.	mm.
EM-8014-UB 0	65	46
EM-8014-UB 1	67	50
EM-8014-UB 2	76	54
EM-8014-UB 3	79	56
EM-8014-UB 4	79	61




For partially toothed upper jaws

	A	B
	mm.	mm.
EM-8014-OP 1	62	50
EM-8014-OP 2	70	56



For partially toothed lower jaws



	A	B
	mm.	mm.
EM-8014-UP 1	70	55
EM-8014-UP 2	78	59



For narrow upper jaws

	A	B
	mm.	mm.
EM-8014-BS 2	70	54


9. Abdruckklöffel · Impression Trays

Endodontic Perforated Without Rimlock



For edentulous upper jaws

	A mm.	B mm.
EM-8014-OU 1	64	46
EM-8014-OU 2	68	53
EM-8014-OU 3	79	61



For edentulous lower jaws

	A mm.	B mm.
EM-8014-UU 1	70	48
EM-8014-UU 2	70	53
EM-8014-UU 3	72	57




9. Abdruckklöffel • Impression Trays

Regular Solid Leaf Handle With Rimlock



Internal high grip satin finish



*Increased diameter (2.1 mm.)
retention rim to provide extra grip.*



EM-8015S-U1

EM-8015S-U2

EM-8015S-U3



EM-8015S-U5

EM-8015S-U5



		A mm.	B mm.
EM-8015S-U1	XL	81	63
EM-8015S-U2	L	74	60
EM-8015S-U3	M	71	55
EM-8015S-U4	S	67	53
EM-8015S-U5	XS	65	51

		A mm.	B mm.
EM-8014S-L1	XL	80	63
EM-8014S-L2	L	78	59
EM-8014S-L3	M	76	57
EM-8014S-L4	S	72	57
EM-8014S-L5	XS	69	50



EM-8014S-L1

EM-8014S-L2

EM-8014S-L3

EM-8014S-L4

EM-8014S-L5




9. Abdrucklöffel · Impression Trays

Regular Perforated With Rimlock



*Increased diameter (2,1 mm.)
retention rim to provide extra grip*



Internal high grip satin finish



EM-8016S-U6


EM-8016S-U5

EM-8016S-U4

EM-8016S-U3



EM-8016S-U2

EM-8016S-U1



		A mm.	B mm.
EM-8016S-U1	XXS	61	47
EM-8016S-U2	XS	64	51
EM-8016S-U3	S	66	53
EM-8016S-U4	M	69	57
EM-8016S-U5	L	71	60
EM-8016S-U6	XL	73	62

		A mm.	B mm.
EM-8016S-L1	XXS	55	47
EM-8016S-L2	XS	60	49
EM-8016S-L3	S	72	54
EM-8016S-L4	M	74	57
EM-8016S-L5	L	76	59
EM-8016S-L6	XL	78	61



XL

L

M

S

XS

XXS



EM-8016S-L6

EM-8016S-L5

EM-8016S-L4

EM-8016S-L3

EM-8016S-L2

EM-8016S-L1




9. Abdrucklöffel · Impression Trays

Regular Solid With Rimlock



Internal high grip satin finish



Increased diameter (2,1 mm.)
retention rim to provide extra grip





9. Abdruckklöffel · Impression Trays

Regular Solid With Rimlock

EM-8017S-U6


EM-8017S-U5

EM-8017S-U4

EM-8017S-U3




EM-8017S-U2

EM-8017S-U1



		A mm.	B mm.
EM-8017S-U1	XXS	61	47
EM-8017S-U2	XS	64	51
EM-8017S-U3	S	66	53
EM-8017S-U4	M	69	57
EM-8017S-U5	L	71	60
EM-8017S-U6	XL	73	62

		A mm.	B mm.
EM-8017S-L1	XXS	55	47
EM-8017S-L2	XS	60	49
EM-8017S-L3	S	72	54
EM-8017S-L4	M	75	57
EM-8017S-L5	L	76	59
EM-8017S-L6	XL	78	61



EM-8017S-L6

EM-8017S-L5

EM-8017S-L4

EM-8017S-L3

EM-8017S-L2

EM-8017S-L1



9. Abdrucklöffel · Impression Trays

Non-Welding Micro Sandblast Trays

	Shape Upper	Width	Length
EM-8018-U2	S	66	53
EM-8018-U3	m	69	57
EM-8018-U4	L	72	60
EM-8017-U5	XS	62	50
EM-8017-U6	XL	75	63

Specialty

High Quality Stainless Steel
Adjustable Front Wale
Sterilizing Autoclave at 134 C
No Welding
Rust-proof
Micro Sand Blast Finish



	Shape Lower	Width	Length
EM-8018-U2	S	72	53
EM-8018-U3	m	75	57
EM-8018-U4	L	78	63
EM-8018-U5	XS	60	49
EM-8018-U6	XL	81	63



Regulars maxilla

EM-4970-110 **EM-4970-111** **EM-4970-112**
 SUP B0 SUP B1 SUP B2

EM-4970-113 **EM-4970-114**
 SUP B3 SUP B4



Edentulous maxilla

EM-4970-131 **EM-4970-132** **EM-4970-133**
 SUP U1 SUP U2 SUP U3



Partials, depressed centers

EM-4970-121 **EM-4970-122** **EM-4970-123**
 SUP P1 SUP P2 SUP P3



Functional impressions

EM-4970-141 **EM-4970-142** **EM-4970-143**
 SUP F1 SUP F2 SUP F3



Regulars mandible

EM-4970-210 **EM-4970-211** **EM-4970-212**
 INF B0 INF B1 INF B2

EM-4970-213 **EM-4970-214**
 INF B3 INF B4



Partials, depressed centers

EM-4970-221 **EM-4970-222** **EM-4970-223**
 INF P1 INF P2 INF P3



Edentulous

EM-4970-231 **EM-4970-232** **EM-4970-233**
 INF U1 INF U2 INF U3



Regulars maxilla

- EM-4972-110** SUP B0 **EM-4972-111** SUP B1 SUP B2
EM-4972-113 SUP B3 **EM-4972-114** SUP B4



Edentulous maxilla

- EM-4972-131** SUP U1 **EM-4972-132** SUP U2 **EM-4972-133** SUP U3



Partials, depressed centers

- EM-4972-121** SUP P1 **EM-4972-122** SUP P2 **EM-4972-123** SUP P3



Functional impressions

- EM-4972-141** SUP F1 **EM-4972-142** SUP F2 **EM-4972-143** SUP F3



Bezahnte Unterkiefer (BU) · Regulars mandible

EM-4972-210 **EM-4972-211** **EM-4972-212**
 INF B0 INF B1 INF B2

EM-4972-213 **EM-4972-214**
 INF B3 INF B4



Partiell bezahnte Kiefer (PU) · Partials, depressed centers

EM-4972-221 **EM-4972-222** **EM-4972-223**
 INF P1 INF P2 INF P3



Unbezahnte Unterkiefer (UU) · Edentulous

EM-4972-231 **EM-4972-232** **EM-4972-233**
 INF U1 INF U2 INF U3



EM-4976-2
FIG. L2
A = 50 mm

EM-4977-2
FIG. R2
A = 50 mm

EM-4978-1
FIG. L1
A = 46 mm

EM-4979-1
FIG. R1
A = 46 mm

EM-4976-3
FIG. L3
A = 55 mm



EM-4977-3
FIG. R3
A = 55 mm

EM-4978-2
FIG. L2
A = 50 mm

EM-4979-2
FIG. R2
A = 50 mm

EM-4978-3
FIG. L3
A = 55 mm

EM-4979-3
FIG. R3
A = 55 mm



EM-4980
STOLLEY
with bolt

EM-4981
STOLLEY
with screw

EM-4982
STOLLEY
with bolt

EM-4983
STOLLEY
with screw

